# Blepha<sup>Plus</sup><sup>®</sup>

**Preservative Free**

**28**

Sterile single-use gauzes for adults and children's pericocular hygiene

## INDICATIONS

Blepha Plus<sup>®</sup> is a preservative free medical device for adults and children particularly adapt for pericocular cleansing, in case of pre and post ophthalmic surgery as well. It is effective in case of inflammations, infections, allergies, ocular secretions, infective or allergic manifestations due to pathologies and for cleansing the eye area.

## PROPERTIES

Hyaluronic acid sodium salt is one of the fundamental substance making up the extracellular matrix of most of tissues. Thanks to its water retention capacity and its filmogen power, it exerts an hydrating and protective activity.

The functional substances cleanse carefully the ocular and periocular area making Blepha Plus<sup>®</sup> a very good hygienic device which can be used also during ocular therapies prescribed by one own's doctor.

The peculiar cleansing, hydrating, protective and smoothing activity is guaranteed by Blepha Plus<sup>®</sup>'s composition resulting to be an extremely gentle product, in particular for more sensitive skin and area.

## USE AND ADMINISTRATION

Wash your hands properly. Open the sachet and extract the gauze. Cleanse carefully the

pericocular area with the gauze in order to remove some ocular secretions or any other residuals.

## INGREDIENTS

100% cotton gauze embedded with 3 ml of the following solution: aqua, PEG/PPG – 20/15 dimethicone, chamomilla recutita flower extract, propylene glycol, euphrasia officinalis extract, sodium hyaluronate, ammonium glycyrrhizate, caprylyl/capryl glucoside, sodium cocoyl wheat amino acids, ethylhexylglycerin phenethyl alcohol, sodium chloride, sodium bicarbonate.

## WARNINGS

- Single use.
- Keep out of the reach of children.
- The sterility of the sachet is guaranteed only in case of undamaged package.
- Store in a cool dry place.
- External use.
- Suspend the treatment in case of clear irritating events or hyper-sensitivity to the product.
- The action and the safety of the product are guaranteed exclusively in case of undamaged package.
- Read carefully the package leaflet.

**NTC SRL**  
Via Luigi Razza 3  
20124 Milano, Italy

**DISTRIBUTOR**  
RockMed Pharma, P.O. Box 121,  
5688 ZJ Oirschot, The Netherlands

[WWW.ROCKMEDPHARMA.NL](http://WWW.ROCKMEDPHARMA.NL)

**ROCKMED**pharma  
innovation education cooperation

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