**Uriflex®**

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USER MANUEL

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MANUEL

INSTRUCCIONES

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NL

EN

DE

FR

ES

IT

**Mickey**

**DE BETERE PLASWEKKERTRAINING!**

**Uriflex®**

Mennonietenbuurt 29 • 1427 AX • Amstelhoek  
Postbus 525 1420 CA • Uithoorn • Tel: 0297 - 760 001  
[www.uriflex.com](http://www.uriflex.com) • [info@uriflex.com](mailto:info@uriflex.com)

CE

## URIFLEX PLASWEKKER "MICKEY"

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## URIFLEX MICKEY BEDWETTING ALARM

Congratulations on your purchase of a Mickey bedwetting alarm and thank you for the confidence you show in this Uriflex bedwetting alarm training!

Since 1954, Uriflex has been specialised in providing bedwetting alarm trainings to prevent bedwetting. In addition to offering good quality wetting alarms Uriflex stands out because of its personal and professional guidance. With our approach approximately 90% of the users will stop bedwetting within 3 months!

It's crucial that the user is motivated to overcome bedwetting.

The bedwetting alarm is not a miracle solution, but offers a training opportunity to also get complete control over the sphincter of the bladder while you sleep. This training can be intense the first few weeks and requires above all teamwork between the supervisor and user, and - if you wish – expert guidance from us. Please contact your local Uriflex branch for more information.

We wish you much success using  
the bedwetting alarm training!

# CONTENT BEDWETTING ALARM SET

- 1 Receiver
- 2 Transmitter
- 3 Sensor briefs \*
- 4 Adapter

- 5 User manual
- 6 Scorecard
- 7 Washing instructions

\* These must be purchased separately.

## HOW TO USE THE MICKEY BEDWETTING ALARM

# HOW TO USE THE MICKEY BEDWETTING ALARM

## To switch on

Plug the connector into both the alarm (see photo) and into a plug socket. On top of the alarm on the left a small light lights up continuously.

## Transmitter

The transmitter includes a 3V-battery. This battery works for at least 2 months in situations of normal use. Before use, all the transmitter press buttons should be snapped onto the front of the sensor briefs (see photo).

- Always remove the transmitter from the sensor briefs when the bedwetting alarm isn't being used. This will spare the battery!
- Replace the battery of the transmitter when the low battery indicator on the bottom of the receiver lights up.

You can always check the proper functioning of the transmitter by connecting the two push buttons on the transmitter using a metal object - for example a pair of scissors. But make sure that the receiver is on. If the alarm doesn't start and the error indicator on the bottom of the receiver doesn't light up, the battery should be replaced.

## Replacing the transmitter battery

Open the transmitter using a flat screwdriver to carefully disconnect the top and bottom part of the transmitter box. Place the screwdriver in the middle of one of the two short sides. Gently remove the present battery from underneath the clamp. Insert a new battery the same way and gently press the clamp. Make sure that the rubber O-rings are placed correctly (thin below and thick above) and shut the transmitter.

## Battery type

The transmitter is powered by a 3V battery type CR2032. We prefer the brand Varta. Most important is that the bottom of the battery is smooth and has as little structure as possible.

## **Connecting the sensor briefs to the bedwetting alarm**

The Mickey bedwetting alarm is wireless, so when the receiver is on and the transmitter is clicked onto the sensor briefs, the bedwetting alarm is ready for use.

## **Select the desired alarm signal**

The Mickey offers 5 different alarms. At the back of the receiver you will find a button with the text "signal". Pressing this button, you will hear an alarm sound. By pressing the button again, the signal changes automatically. This way you can switch between up to 5 different signals. To select the desired alarm release the button when you hear the signal of your choice.

## **Setting the volume**

When pressing the volume button on the bottom of the receiver, you will hear the volume as it is set on the bedwetting alarm. By keeping the button pressed for a while the volume changes automatically. To select the desired volume release the button when you hear the volume of your choice.

## **Low battery signal**

On the back of the receiver you will find a light that says "low battery". If this (orange) light lights up, the battery of the transmitter must be replaced. For instructions on how to open the transmitter, see previous page.

## **Error signal**

When a red "error" light lights up on the back of the receiver this is an indication that the transmitter and the receiver are not connected properly. When wetting the bed the alarm will not go off.

### **In this case, please follow the following steps:**

1. Remove the battery from the transmitter
2. Unplug the receiver from the plug socket
3. Wait 1 minute
4. Plug the receiver back into the plug socket
5. Place the battery back into the transmitter

If this doesn't help, you should repeat the steps 1 to 4. For step 5 please use a new battery. Please contact our local office (see [www.uriflex.com](http://www.uriflex.com)) if the error light still lights up despite the new battery.

### **Bedwetting alarm ready for nightly use**

At bedtime, the Mickey bedwetting alarm is switched on, the Uriflex sensor briefs are put on and the transmitter is attached to them. Place the receiver within the range of a few meters at a location in the bedroom that you like best. The only condition is that the user really needs to get out of bed to turn off the alarm! The user must wake entirely and should not be able to switch off the alarm subconsciously.

### **Then the alarm goes off...**

As soon as the first drops of urine reach the briefs, the alarm will sound immediately. The user awakes up and in a reflex, the sphincter of the bladder closes. The user gets out of bed, first unplugs the transmitter from the sensor briefs and subsequently switches off the alarm by pressing the red button on top of the receiver. Now, the user finishes urinating on the toilet. After putting on a pair of clean sensor briefs and snapping the transmitter back on, the user can drift back off to sleep.

### **Important: All of this must happen in a conscious state!**

It is crucial that the user is awake and fully aware of all the actions taking place, from waking up and urinating on the toilet to putting on a clean pair of sensor briefs. It may be desirable that a parent or guardian temporarily sleeps in the same room as the user, to make to sure the user immediately responds to the alarm.

Immediate responding to the alarm is definitely the key to success. Only then will the necessary reflex of waking up be encouraged and developed.

The user must wake up because of the alarm. The alarm only sounds when the first drops of urine reach the sensor briefs. This is the critical moment the user is training to get ahead of. If, for example, you wake up the user at any time during the night to urinate, he/she won't develop an automatic reflex. By using the Uriflex bedwetting alarm every night, the user will start waking up on time by him/herself after already a few weeks.

### **Completing the scorecard together every morning**

Every morning a sticker is placed on the supplied scorecard and you write down how many points have been earned that night. Or you write down how well the night went. This of course has a very stimulating effect,

because the user wants to earn as many points as soon as possible. This can be extra effective if you link a small personal reward to these accomplishments.

If you send us your completed scorecard, we can monitor the progress of the training. And if you wish, we can always provide you with personalised advice on the basis of the progress.

The training course is successful if the user stays dry all night for three consecutive weeks or if he/she wakes up without help from the alarm to go to the toilet.

### **Personal help desk**

We can personally guide you through the course of the bedwetting alarm training. See [www.uriflex.com](http://www.uriflex.com) for details of your local branch.

### **Practical tips**

- For very restless sleepers it's recommended put on their own underwear on top of the Uriflex sensor briefs to avoid the transmitter coming off during the night.
- If the user perspires a lot, it's advisable to put on normal underwear under the Uriflex sensor briefs. The alarm will then not sound after the first drops of urine, but it will prevent unnecessary alarms.

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