

RESMED

Simple to use,
easy to fit

www.resmedmiragefx.com

RESMED

Simple to use,
easy to fit

User Information

Easy to fit

Simple to use



To get the most out of your therapy, pair your Mirage FX or Mirage FX for Her with ResMed's new S9™ Series of therapy devices. Visit www.resmed.com for more information.

ResMed Ltd.
1 Elizabeth Macarthur Drive
Bella Vista NSW 2153 Australia.

Distributed by: ResMed Corp 9001 Spectrum Center Boulevard San Diego CA 92123 USA.
EC/REP ResMed (UK) Ltd 96 Milton Park Abingdon Oxfordshire OX14 4RY UK. See www.resmed.com for other ResMed locations worldwide. S9, Mirage, SoftEdge and Spring Air are trademarks of ResMed Ltd. S9, Mirage and Spring Air are registered in U.S. Patent and Trademark Office. © 2011 ResMed Ltd Specifications may change without notice. 1014962/1 2011-07

Mirage™ FX
NASAL MASK


Mirage™ FX for Her
NASAL MASK

Mirage™ FX
NASAL MASK

Mirage™ FX for Her
NASAL MASK

Simple to use, **easy to fit**

Beyond the simplicity of its minimalist design is intelligent technology working hard to provide the best therapy and user experience possible.



It's all **about you**

Mirage FX is so streamlined and unobtrusive, you'll love the feeling of openness and clear line of sight ... and being one of the lightest nasal masks around, Mirage FX won't weigh you down.

Made up of only four parts, including headgear, Mirage FX is easy to clean, easy to fit and easy to adjust, without compromising quality, comfort or performance.

Mirage™ FX
NASAL MASK

www.resmedmiragefx.com

Exciting features, **real benefits**



- 1 Form-fitting SoftEdge™ headgear**
Springs to life, ready to fit and wear; is kind to your face, reducing marks and discomfort
- 2 Flexible, winged forehead support**
Allows a dynamic range of fit and adjustment
- 3 Lightweight, durable frame**
Strong enough to stabilise the mask on your face, soft enough to accommodate your movement during sleep
- 4 Dual-wall, Spring Air™ cushion**
Provides even pressure distribution, with extra softness at the nasal bridge, reducing common discomfort in this sensitive area
- 5 Diffuse vent design**
Disperses air gently and quietly—at only 26 dbA*—away from your bed partner, for a good night's sleep for all
- 6 Headgear loops**
Make it easy to attach and detach the headgear, while retaining optimal headgear settings
- 7 Squeeze-tab elbow**
Makes it easy to attach and detach the elbow from the frame

Something for **everyone**

The Mirage FX range of masks means there is a solution for everybody.

Mirage FX: the standard choice

Using thousands of anthropometric facial scans, the standard Mirage FX size was designed to cater for a remarkable 9 out of 10 people.

Mirage FX Wide: a broad fit range

To ensure that everyone can achieve maximum fit and comfort, Mirage FX Wide caters for those with wider nasal features.

Mirage FX for Her: just for you

Designed exclusively for women, Mirage FX for Her provides all the same features and benefits as Mirage FX, but with a feminine touch.

The headgear is smaller with a snug crown strap and hair management options. Whether you have long or short hair, wear it tied or untied, the custom-made headgear back-strap easily accommodates your personal style... and it's finished in a soft and stylish colour scheme.



Mirage™ FX for Her
NASAL MASK

* Testing in accordance with ISO 4871.