

Breeze Nasal

1 . Silicone Cushion

Its soft silicone material provides optimum comfort. It easily snaps onto the rigid shell. 3 sizes of cushions are available to fit all morphologies (S, M, L).

2 . Frame

Light and durable, it allows for quick assembly and disassembly of the silicone cushion and the frame assuring better hygiene.

3 . Forehead support tab

Made of flexible silicon, it provides a comfortable mask fit while relieving pressure points on the nasal bridge.

4 . Rotating Elbow

360° rotating elbows provide greater latitude for the patient to move naturally throughout the night.

5 . Noiseless

Designed for quiet exhalation, ensuring a peaceful sleep for the patient and the patients bed partner.

6 . Headgear

Rapid Velcro® fitting system and the headgear hook allow for an optimal fit of the mask, while keeping the headgear settings as chosen.



Available in 3 sizes