

# Breeze Comfort Range

The Breeze Comfort Range offers an innovative concept.

Single hand setting of the forehead support ensures simple mask adaptation for all patients and the quick coupler system's headgear clips allow for easy fitting and removal of the mask.


Set it once and forget it !

## Facial Comfort



Available in 3 sizes

## Nasal Comfort



Available in 3 sizes

### 1 . Silicone Cushion

Deep silicone bubble offers maximum comfort for patients with a prominent nose.  
The surface of the mask in contact with the skin is supple.

### 2 . Frame

Light and durable, it allows for quick assembly and disassembly of the silicone cushion and the frame assuring better hygiene. .

### 3 . Adjustable Forehead support

The intuitive adjustment of the frontal support allows for optimal tensioning of the headgear .

### 4 . Rotating Elbow

360° rotating elbows provide greater latitude for the patient to move naturally throughout the night.

### 5 . Noiseless

Designed for quiet exhalation, ensuring a peaceful sleep for the patient and the patients bed partner.

### 6 . Headgear

Rapid Velcro® fitting system and the quick attach clips allow for an optimal fit of the mask, while keeping the headgear settings as chosen.