


Breeze Pillows Mask



1 . Nasal pillows

Packed with three pillow sizes to suit most morphology.
Anatomical forms allow for ideal positioning.

2 . Silent

The vents are designed to diffuse the air in silence, ensuring a peaceful sleep.

3 . Rotating elbow

360° rotating elbows provide greater latitude for the patient to move naturally throughout the night.

4 . Headgear

Easy to set up and remove.
One size adjustable easily.
A strap for positioning the silicon tube on the top of the headgear.

5 . Silicone tube

Lightweight and flexible, it reduces pull on the mask enhancing stability and comfort.

With minimal contact on the face, the Breeze Pillows Mask, offers a real sense of freedom .

The Breeze Pillows Mask has all the advantages to facilitate good night's sleep, allowing greater freedom of movement, while ensuring stability of the mask through the night.