

# Breeze

## Nasal



*Available in 3 sizes*

### 1 . Silicone Cushion

Deep silicone bubble offers maximum comfort for patients with a prominent nose.

The surface of the mask in contact with the skin is supple .

### 2 . Frame

Light and durable, it allows for quick assembly and disassembly of the silicone cushion and the frame assuring better hygiene .

### 3 . Forehead support tabs

For optimum adjustment, 3 sizes of interchangeable forehead support tabs are supplied with the mask.

### 4 . Rotating Elbow

360° rotating elbows provide greater latitude for the patient to move naturally throughout the night.

### 5 . Noiseless

Designed for quiet exhalation, ensuring a peaceful sleep for the patient and the patients bed partner.

### 6 . Headgear

Rapid Velcro® fitting system and the headgear hook allow for an optimal fit of the mask, while keeping the headgear settings as chosen.



## Facial

*Available in 3 sizes*