

Dental protection against wear caused by bruxism **Thermo-formable, flexible and comfortable**

Bruxism is movement of the teeth caused by either clenching or grinding. It can be caused by stress and can lead to wear of the teeth, and even fractures and joint problems.

Principle of operation

Thermo-formable, flexible dental splint. Protects the teeth from wear caused by bruxism, grinding and clenching of the teeth.



Indications

- Protects teeth against erosion due to bruxism, grinding and clenching of the teeth.

Contraindications

- Severe periodontal disorders.
- Fewer than 8 teeth on the lower arch.
- If the teeth are unstable, ask your dentist for advice.
- Users under 18 years of age.

Contents

Ref. S-04-PRO

- One flexible thermoformable dental splint for the lower jaw.
- One toothprint plate.
- One floating thermometer.
- One ventilated storage box.

Patented medical device
CE

