

Positioning pillow Anti-snoring

According to scientific studies, snoring for 55% of snorers is due to a positional problem. They snore more sleeping on the back than on the side. POSIFORM has a unique anti-snoring design, which ensures maximum effectiveness.

Principle of operation


This pillow has a form which stimulates sleeping on the side with the head inclined slightly downward. In this position, the tongue is pushed slightly forward by gravity disengaging the pharynx and ensuring unobstructed breathing.

Indications

- Positional snoring.
- Light and positional apnoea; requires authorisation from the physician.

Contraindications

- Articular neck pains.
- Very light sleep with frequent movements.



Effectiveness

Very effective in the case of positional snoring and stable sleep.

Contents

Ref. S-PSF

- 1 memory foam contour pillow.
- 1 with interior cotton cover and bamboo cover.

