

# VICRYL\* Plus

(Polyglactin 910) Braided Antibacterial

## SUTURE

<b>Description:</b>	<b>Braided</b>	For easy handling and secure knot tying.
	<b>Coated</b>	For smooth passage through tissue and easy knot tie down.
	<b>Synthetic</b>	For minimal tissue reaction.
	<b>Absorbable</b>	By hydrolysis. Essentially complete in 56-70 days.

<b>Tensile strength:</b>	<b>Post implantation</b>	<b>Approximate % original strength remaining</b>
	14 days	75%
	21 days	50%
	28 days	25%

**Color:** Violet or undyed.

**Range:** 5-0 to 2 (USP). Supplied as needled sutures and ligatures.

**Indications:** VICRYL\* Plus sutures are intended for use in soft tissue approximation and/or ligation, including use in microsurgery for vessels less than 2 mm in diameter. The safety and effectiveness of VICRYL\* Plus sutures in cardiovascular tissue, ophthalmic surgery and neurological tissue have not been established. Typical areas of use include:


- Fascia closure - due to the 28-day claim
- Subcutaneous fat
- Joint capsule
- Uterus





### Protection against bacterial colonization


In vitro studies demonstrate, that VICRYL\* Plus suture has a zone of inhibition that is effective against the pathogens that most frequently cause surgical site infection, including *Staphylococcus aureus*, methicillin resistant *Staphylococcus aureus* (MRSA), *Staphylococcus epidermidis* and methicillin resistant *Staphylococcus epidermidis* (MRSE).


Full details in the Instruction for Use included in every package.



3/8 Circle Reverse Cutting PRIME* Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 PS-1 24 mm	3-0	○	45		36	VCP683H
	3-0	○	45	MULTIPASS*	36	MPVCP683H
	2-0	○	45		36	VCP684H

3/8 Circle Reverse Cutting Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 FS-2 19 mm	5-0	○	45		36	VCP291H
	5-0	●	45		36	VCP391H
	4-0	○	45		36	VCP292ZH
	4-0	●	45		12	VCP392ZG
	4-0	●	45		36	VCP392ZH
	4-0	○	70		36	VCP422H
	4-0	●	70		36	VCP397H
	3-0	○	45		36	VCP293H
	3-0	●	45		36	VCP393H
	3-0	○	70		36	VCP423H
	3-0	●	70		36	VCP398H
 FS-2S 19 mm	4-0	○	45	Slimmed	36	VCP292H
	4-0	●	45	Slimmed	36	VCP392H
 FS-1 24 mm	4-0	○	70		36	VCP441H
	4-0	●	70		36	VCP451H
	3-0	○	70		36	VCP442H
	3-0	●	70		36	VCP452H
	2-0	○	70		36	VCP443H
	2-0	●	70		36	VCP453H
 FSL 30 mm	3-0	●	70		36	VCP585H
	2-0	○	70		36	VCP589H
	2-0	●	70		36	VCP586H
	0	●	70		36	VCP587H








3/8 Circle Conventional Cutting PRIME* Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 PC-3 16 mm	5-0	○	45		36	VCP9505H
	4-0	○	45		36	VCP9506H

Straight Cutting Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 KS 60 mm	3-0	○	70		36	VCP663H





# VICRYL\* Plus

(Polyglactin 910) Braided Antibacterial

## 1/2 Circle Reverse Cutting Needles (Heavier Wire Diameter)

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 J-1	17 mm	4-0	●	75		36	VCPJV457H
 X-1	22 mm	3-0	●	45		36	VCP460H
		3-0	○	70		36	VCP458H
		2-0	○	70		36	VCP459H
 OS-4	22 mm	0	○	70	Heavy Needle	36	VCP694H
		1	○	70	Heavy Needle	36	VCP695H
		2	●	70	Heavy Needle	36	VCP9287H
 CP-2	26 mm	2-0	●	70		36	VCP969H
		0	●	70		36	VCP974H
 CP-1	36 mm	2-0	●	70		36	VCP466H
		2-0	○	90		36	VCP1472H
		0	●	70		36	VCP467H
		0	●	90		36	VCP473H
		1	●	70		36	VCP468H
		1	●	90		36	VCP474H
		2	●	90		36	VCP482H
 CP	40 mm	2-0	●	90		36	VCPJV484H
		0	●	70		36	VCP479H
		0	●	90		36	VCP485H
		1	●	70		36	VCP480H
		1	●	90		36	VCP486H
 CPX	48 mm	0	●	70		36	VCP9295H
		0	●	90		36	VCPV1057H
		1	●	90		36	VCP1058H
		2	●	90		36	VCP1059H

## 3/8 Circle Reverse Cutting PRIME\* Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 P-3	13 mm	5-0	○	45	MULTIPASS*	36	MPVCP493H
		5-0	○	45		36	VCP493H
		4-0	○	45	MULTIPASS*	36	MPVCP494H
		4-0	○	45		36	VCP494H
 PS-3	16 mm	5-0	○	45	MULTIPASS*	36	MPVCP500H
		5-0	○	45		36	VCP500H
 PS-2	19 mm	4-0	○	45	MULTIPASS*	36	MPVCP496H
		4-0	○	45		36	VCP496ZH
		3-0	○	45	MULTIPASS*	36	MPVCP497H
		3-0	○	45		36	VCP497H
 PS-2 S	19 mm	4-0	○	45	Slimmed	36	VCP496H